DK	Blodtryksmålør.....	2
SE	Blodtrycksmätare.....	9
NO	Blodtrykksmålør.....	16
FI	Verenpainemittari.....	23
UK	Blood pressure monitor.....	30
DE	Blutdruckmessgerät.....	37
PL	Cięśnieniomierz.....	44

INTRODUKTION

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INTRODUCTION

To get the best out of your new appliance, please read this user guide carefully before using it for the first time. Take particular note of the safety precautions. We also recommend that you keep the instructions for future reference.

IMPORTANT SAFETY PRECAUTIONS

- Incorrect use of this appliance may cause personal injury and damage.
- Use for its intended purpose only. The manufacturer is not responsible for any injury or damage resulting from incorrect use or handling (see also Guarantee Terms).
- For domestic use only. Not for outdoor or commercial use.
- This appliance must not be used to make medical diagnoses. It is suitable only for reference. Consult a doctor if your blood pressure readings seem abnormal.
- Always follow the procedure specified in these user instructions when taking readings.
- Do not stretch, twist or bend the cuff and hose.
- Do not use on babies or others who cannot communicate if the cuff is too tight, etc.
- Do not immerse in water or any other liquid.
- Do not use in the vicinity of mobile phones or microwave ovens as these emit electro-magnetic radiation that can affect the appliance.
- Do not expose to direct sunlight, high temperatures, humidity, noise or corrosive substances.
- Do not use if not functioning properly, or if damaged by a heavy impact or similar.


- If the appliance, hose or cuff have been damaged, have them inspected and if necessary repaired by an authorised repair engineer. Never try to repair the appliance yourself. Please contact the store where you bought the appliance for repairs under guarantee. Unauthorised repairs or modifications to the appliance will invalidate the guarantee.

Batteries

- Remove the batteries if the appliance is not to be used for a prolonged period.
- Change the batteries immediately once they are dead.
- Never use different types of batteries, or use a new and an old battery at the same time.
- Ensure the batteries are correctly polarised when installing them.


KEY TO MAIN COMPONENTS OF THE APPLIANCE

1. Velcro strap
2. Cuff
3. Ø symbol
4. Air hose
5. Air hose hole
6. Display
7. Control buttons
8. Hose connector
9. <-OK mark
10. OK-> mark
11. INDEX mark

Display and control buttons

- a. Systolic blood pressure
- b. Diastolic blood pressure
- c. Pulse
- d. On/off button
- e. M2 button (memory 2)
- f. M1 button (memory 1)
- g. Battery symbol
- h. Pulse symbol

ABOUT BLOOD PRESSURE

What is blood pressure?

When the heart pumps blood around the body, the blood exerts pressure on the walls of the blood vessels. This is known as blood pressure.

Blood pressure is different in different parts of the body and depends on the pulse, how elastic the blood vessels are and how thick the blood is.

The pressure that can be measured when the heart pumps out into the body is known as systolic blood pressure. The pressure that can be measured when blood runs back from the body towards the heart is known as diastolic blood pressure.

The two types of blood pressure are usually given together, with the systolic blood pressure always given first. If blood pressure is given as "120/70" (expressed as "120 over 70"), this means that the systolic blood pressure is 120, and the diastolic 70.

The reading unit for blood pressure is mmHg (millimetres of mercury).

Why does blood pressure vary?

Blood pressure varies throughout the day and over a person's life. Many different factors affect your blood pressure, something you should be aware of when taking your blood pressure. For instance, be aware of the following:

- Blood pressure is higher during the day than at night.
- Blood pressure increases slightly after you have just eaten. A lot of salt in food can also increase blood pressure.
- Light exercise (e.g. walking and jogging) causes the systolic blood pressure to rise slightly, while hard exercise will cause it to rise dramatically. However, it will fall again when the body relaxes. A lack of exercise raises blood pressure permanently.
- Blood pressure can increase if you are excited, stressed or tense.
- Sudden climatic changes can cause blood pressure to rise.

- The systolic pressure increases with age.
- Women have a greater risk of increased blood pressure after the menopause.
- Insufficient sleep, smoking, obesity, alcohol and certain illnesses cause blood pressure to rise.
- Bathing and toilet visits may cause blood pressure to rise.

When is my blood pressure too high or too low?

High blood pressure is also called hypertension. Low blood pressure is also called hypotension. The table below shows general guidelines for when your blood pressure is too high or too low. You should always consult your doctor for exact readings.

Range	Systolic	Diastolic
Very high blood pressure (hypertension)	> 140 mmHg	> 90 mmHg
High normal	130-139 mmHg	85-89 mmHg
Normal	91-129 mmHg	51-84 mmHg
Very low blood pressure (hypotension)	< 90 mmHg	< 50 mmHg

Why should I measure my blood pressure?

Blood pressure rises naturally with age. But high blood pressure can also be the first sign of illness, being overweight or of the arteries beginning to stiffen, which increases the risk of blood clots, among other things. It is therefore a good idea to keep an eye on your blood pressure and consult a doctor if it seems unusually high or low. Write down your own readings, or save them using the memory function, to be able to show them to your doctor.

Important! Do not use to make medical diagnoses. For use only as reference. Consult a doctor if your blood pressure readings seem abnormal. Always follow the procedure specified in these user instructions when taking readings.

PRIOR TO FIRST USE

Inserting the batteries

- Caution! To save readings in the memory, change the batteries within one minute. If changing the batteries takes more than

one minute, the display will show all its symbols for a moment, and will require resetting (all saved readings will be deleted).

- Change the batteries when the battery symbol (g) is shown in the display (6).
- Remove battery cover by pressing the lock tab and lift up.
- Insert four AAA batteries. Check batteries are correctly polarised (see diagram in battery compartment.)
- Replace the cover.

Fitting cuff and air hose

- Connect hose connector (8) to air hose hole (5).
- Check connector is pushed in as far as possible and is secure.

USING THE APPLIANCE

Preparing for a blood pressure reading

- Sit or lie down for 5-10 minutes before the reading.
- Avoid eating, smoking, exercising, bathing or drinking alcohol, coffee or tea before taking the reading.
- Do not go to the toilet immediately before the reading.
- Do not take a blood pressure reading if you are excited, stressed or tense.

- Relax during the reading and do not talk.

Blood pressure reading

- Rest arm on a cushion or table top with palm upwards.
- Place cuff (2) on upper arm, approx. 2-3 cm above the elbow.
 - o The cuff should be approximately level with your heart to obtain an accurate reading.
- Turn the cuff so that the "Ø" (3) mark just under the air hose (4) is on the inside of arm, facing downwards.
- Slip your free hand under the cuff and grasp the Velcro strap (1) at the "INDEX" mark (11).
- Pull strap around cuff, and bind tightly.
- When cuff is securely on arm, the "INDEX" mark should be between "OK->" and "<-OK" (9 and 10). If this is not possible, the cuff is too small for your arm.
- NB:
 - o The cuff should be taut but should not pinch.
 - o The "Ø" mark must be on the inside of the arm, over the major arteries.
 - o The side of the cuff facing your hand must be 2-3 cm from the hand.
- Switch on using the on/off button (d).
- All symbols and figures will flash on the display.
- Pump cuff up to approx. 165 mmHg. The pressure will be constantly shown on the display.
 - o The appliance will hum while the cuff is pumped up and tightens.

- When the cuff is pumped up to approx. 165 mmHg, the figures on the display will drop, and blood pressure reading starts.
 - o The appliance will show the pressure continuously on the display while reading is in progress.
 - o The pulse is also measured, and the pulse symbol (h) is shown on the display.
- When the reading is complete, the systolic pressure (a), diastolic pressure (b) and pulse (c) will be displayed.
- The pulse symbol will flash if the blood pressure reading is too high or too low (see table above). The pulse symbol will remain stable if the reading is within the normal range.
- Results will be shown on the display for approx. 3 minutes and then automatically saved in the appliance memory (M1/memory 1) when the appliance is switched off.
 - o You can also save the result in one of the two memories, by pressing either M1 (f) or M2 (e). You can then switch the appliance off using the on/off button.
- If you need to take several readings, you should wait 5-10 minutes between each reading.
 - o Caution! A single reading will not give a sure indication of your condition. You can only get an idea of this after a number of readings.

Storing blood pressure readings

- The appliance can save 100 readings – 50 in each of the two memories (M1 and M2). The latest reading is always saved as number 01 in the memory. Older readings are moved down the scale as new ones are saved, and deleted when the memory is full.

Displaying saved blood pressure readings

- Press M1 or M2 to show saved readings.
 - o The reading number (e.g. "01") will be displayed for approx. one second. The details for that reading will then be displayed: systolic and diastolic blood pressure and pulse.
 - o A small figure 1 or 2 will be displayed (over the pulse symbol) which indicates which memory is active.
- Press the button again to see older readings.

Deleting saved blood pressure readings

- Remove the batteries and wait approx. 5 minutes.
- When replacing the batteries, the display will show all its symbols for a moment, and will require resetting (all saved readings will be deleted). The memory will now be deleted.

ERROR MESSAGES

The appliance registers any problems and shows the following error messages on the display:

Error	Cause	Solution
Error during pumping	<ol style="list-style-type: none">1. The cuff is not correctly attached.2. The cuff is broken.3. The pump or valves are defective.	<ol style="list-style-type: none">1. Attach the cuff correctly.2. Contact the retailer.3. Contact the retailer.
Disturbances or noise during	<ol style="list-style-type: none">1. The user moves or speaks during the reading.2. The user is too nervous.	<ol style="list-style-type: none">1. The user must sit still and avoid talking during the reading.2. Have a blood pressure reading taken by a doctor.
Error while emptying the	<ol style="list-style-type: none">1. The valves are defective.2. The cuff is not correctly attached.	<ol style="list-style-type: none">1. Contact the retailer.2. Attach the cuff correctly.
Abnormal reading	The user may have a heart problem.	Check result by performing one or more new readings. Talk to your doctor if you consistently have abnormal readings.
Cuff pressure too high.	The air tube is blocked, the pump mechanism is defective or the pressure monitor is defective.	Try to take another reading, or contact the retailer.
Flat batteries (battery symbol on display) or blank	The batteries are almost dead.	Replace the batteries. Check that the batteries are in the right way round.

STORAGE

When the appliance is not in use, store in the accompanying storage box.

- Store in a dry, cool place.
- Remove the batteries if the appliance is not to be used for a prolonged period.
- Do not bend cuff or air hose excessively.
- Do not place heavy objects on the appliance.

CLEANING

When cleaning the appliance, you should pay attention to the following points:

- Clean cuff and appliance with a well-wrung cloth with a little detergent if necessary.
- Do not use strong, corrosive or abrasive cleaning agents for cleaning the appliance, as these may scratch the surface.

INFORMATION ON DISPOSAL AND RECYCLING OF THIS PRODUCT

Please note that this Adexi product is

marked with this symbol:

This means that this product must not be disposed of together with ordinary household waste, as electrical and electronic waste must be disposed of separately.

In accordance with the WEEE directive, every member state must ensure correct collection, recovery, handling and recycling of electrical and electronic waste. Private households in the EU can take used equipment to special recycling stations free of charge. In some member states you can, in certain cases, return the used equipment to the retailer from whom you purchased it, if you are purchasing new equipment. Contact your retailer, distributor or the municipal authorities for further information on what you should do with electrical and electronic waste.

GUARANTEE TERMS

The guarantee does not apply:

- if the above instructions are not followed
- if the appliance has been interfered with
- if the appliance has been mishandled, subjected to rough treatment, or has suffered any other form of damage

Due to the constant development of our products in terms of function and design, we reserve the right to make changes to the product without warning.

QUESTIONS AND ANSWERS

If you have any questions regarding use of the appliance and cannot find the answer in this user guide, please try our website at www.adexi.eu.

Go to the "Consumer Service" menu, click on "Question & Answer" to see the most frequently asked questions.

You can also see contact details if you need to contact us for technical issues, repairs, accessories and spare parts.

IMPORTER

Adexi Group
www.adexi.eu

We cannot be held responsible for any printing errors.

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